SCHOOL OF PHYSICS – LEVEL 200 PHYSICS LABORATORY REPORT SLIP FOR ZCT 293/2 AND ZCT 294/2

Instructions to student: Please make sure you fill in the form completely.										
Instructions to lect	curer: Kindly record the numerical marks in the rub	ric asse	essm	ent for	m, no	ot here				
	PARTICULARS									
Name	: TAN WEI LIANG									
Matric no.	: 22302889									
Group	: W2									
Expt. Code	: 2OS2									
Expt. Title	: DIFFRACTION GRATING SPECTROMETER									
Lecturer in charge	: DR. SITI AZRAH MOHAMAD SAMSURI									
Report due date	: 21/12/2024									
	Lab Report Grade (✓)									
1 🗆	2 □ 3 ☑ 4 □ 5 □ 6 □	A		A-		B+		В		
		В-		C+		C		C-		
		D+		D		D-		F		
	Date Received (Stamp):				Con	nments	:			

REPORT SUBMISSION FORM



Name	:	TAN WEI LIANG
Partner's Name	:	TAN SHUENN PYNG
Group	:	W2
Experiment Code	:	2OS6
Experiment Title	:	DIFFRACTION GRATING SPECTROMETER
Lecturer's/Examiner's Name	:	DR. SITI AZRAH MOHAMAD SAMSURI
Starting Date (1st session)	:	20/11/2024
Ending Date (2nd session)	:	27/11/2024
Submission Date	:	21/12/2024

DECLARATION OF ORIGINALITY

- I, **TAN WEI LIANG 22302889** hereby declare that this laboratory report is my own work. I further declare that:
 - 1. The references/bibliography reflect the sources I have consulted, and
 - 2. I also certify that this report has not previously been submitted for assessment in this or any other units, and that I have not copied in part or whole or otherwise plagiarized the work of other students and/or persons.
 - 3. Sections with no source referrals are my own ideas, arguments, and/or conclusions.

	(1_	
Signature:		Date: 21/12/2024

DIFFRACTION GRATING SPECTROMETER

By

TAN WEI LIANG

21 December 2024

DIFFRACTION GRATING SPECTROMETER

ABSTRACT

The title of this experiment is **DIFFRACTION GRATING SPECTROMETER**. In this experiment, the focus was on analyzing the bright-line emission spectra of Mercury, Helium, Cadmium, Hydrogen, and Sodium vapor light sources using a diffraction grating spectrometer. Light emitted by these vapor sources was diffracted through a grating, and the angles of the resulting spectral lines were measured. The calculated wavelengths corresponded closely to the known values for each element, validating the quantized nature of atomic electronic transitions. Mercury showed characteristic lines at 431.48 nm (purple), 541.97 nm (green), and 572.40 nm (yellow). Helium exhibited spectral lines at 452.09 nm (blue), 501.32 nm (cyan), and 580.58 nm (yellow). Cadmium showed distinct lines at 468.09 nm (blue), 480.85 nm (light blue), and 510.57 nm (green). Hydrogen exhibited lines at 431.88 nm (blue), 462.25 nm (cyan), and 661.52 nm (red). Finally, for Sodium, a prominent line was observed at a calculated wavelength of 580.20 nm (orange). Minor discrepancies between the calculated and theoretical values were attributed to experimental limitations, including apparatus alignment, angular measurement errors, and challenges in observing faint spectral lines.

ACKNOWLEDGEMENTS

First and foremost, I express my sincere appreciation to DR. SITI AZRAH MOHAMAD SAMSURI, our distinguished lecturer and examiner, for the invaluable guidance and unwavering support extended throughout our scientific exploration. I extend my sincere gratitude to my experiment partner, TAN SHUENN PYNG. His invaluable cooperation and dedication throughout both experiments were instrumental to the success of this project. I appreciate his commitment, expertise, and teamwork, which made these scientific endeavours both productive and enjoyable.

CONTENTS

A]	BSTRACT	i
A	CKNOWLEDGEMENTS	ii
LI	IST OF TABLES	v
LI	IST OF FIGURES	vi
1	INTRODUCTION	1
2	THEORY	2
3	METHODOLOGY	4
4	DATA ANALYSIS	6
	4.1 Mercury	8
	4.2 Helium	10
	4.3 Cadmium	11
	4.4 Hydrogen	12
	4.5 Sodium	
5	DISCUSSION	14
6	CONCLUSION	16
R	EFERENCES	17
A]	PPENDICES	19
	Data	19
	Mercury	23
	Mercury percentage difference	25
	Helium	27
	Helium percentage difference	29
	Cadmium	31
	Cadmium percentage difference	33
	Hydrogen	34
	Hydrogen percentage difference	36

Sodium	37
Sodium percentage difference	39

LIST OF TABLES

1	Spectral Line Data for Various Elements with θ in Radians	6
2	Spectral Line Data for Mercury (Hg).	8
3	Mercury spectrum data with calculated and real wavelengths, and their percent-	
	age differences.[4][5]	9
4	Spectral Line Data for Helium (He)	10
5	Helium spectrum data with calculated and real wavelengths, and their percent-	
	age differences.[6][7]	10
6	Spectral Line Data for Cadmium (Cd).	11
7	Cadmium spectrum data with calculated and real wavelengths, and their per-	
	centage differences.[8][9]	11
8	Spectral Line Data for Hydrogen (<i>H</i> ₂).	12
9	Hydrogen spectrum data with calculated and real wavelengths, and their per-	
	centage differences.[10][11]	12
10	Spectral Line Data for Sodium (Na).	13
11	Sodium spectrum data with calculated and real wavelengths, and their percent-	
	age differences.[12][13]	13

LIST OF FIGURES

1	Ray diagram for first order diffraction pattern	2
2	Scattering diagrams showing the interaction of light waves with a diffraction	
	grating.	3
3	Experimental setup for the diffraction grating spectrometer. The setup includes	
	the spectrometer (center), a light source (right) emitting the spectrum, and a	
	telescope for observing and measuring diffraction angles. The diffraction grat-	
	ing is positioned between the collimator and the telescope to disperse the light	
	into spectral lines.	4

INTRODUCTION

The purpose of this experiment is to determine the wavelengths of the colors in the emission spectra of sodium, mercury, helium, cadmium, and hydrogen vapor. An emission spectrum is produced when light emitted from an excited gas passes through a diffraction grating. Unlike a continuous spectrum produced by a solid incandescent source, the light emitted from an electric discharge in a rarefied gas of a single element consists of discrete wavelengths, forming a bright-line spectrum. Each line in the emission spectrum corresponds to a unique electronic transition within the atoms of the element, making the pattern of colors characteristic of that element. The light is separated into its spectral components by passing it through a narrow slit and then dispersing it with a diffraction grating. These separated colors, visible as bright lines, provide insight into the atomic structure and energy levels of the element. The focus of this experiment is to analyze and measure the angular positions of these bright lines to calculate their corresponding wavelengths, which will be compared to the known emission lines of sodium, mercury, helium, cadmium, and hydrogen.

THEORY

An incandescent source, such as a hot solid metal filament, produces a continuous spectrum of wavelengths. In contrast, light emitted by an electric discharge in a rarefied gas of a single element contains a limited number of discrete wavelengths, known as an emission or "bright line" spectrum. The pattern of colors in an emission spectrum is characteristic of the element. These individual colors appear in the form of bright lines because the light, which is separated into the spectrum, usually passes through a narrow slit illuminated by the light source [1].




Figure 1: Ray diagram for first order diffraction pattern

A diffraction grating is a transparent material on which a large number of equally spaced parallel lines have been ruled. The distance between these lines is called the grating line spacing, denoted by d. When light strikes the grating, it is diffracted by these parallel lines. The diffracted light passes through the grating at all angles relative to the original light path. If the diffracted light rays from adjacent lines are in phase, an image of the light source is formed. This happens when the difference in path length between the light rays equals an integral multiple of the wavelength of the light.

The diffraction angle θ , at which light rays are diffracted by the grating, depends on the grating line spacing d and the wavelength λ of the light. The relationship is described by the diffraction equation:

$$\lambda = d \sin \theta$$

This equation allows for the determination of the wavelength λ when the diffraction angle

 θ and the grating spacing d are known. In the **Figure 1**, the path difference for Ray A is one wavelength longer than the path difference for Ray B, illustrating the principle of constructive interference at the diffraction grating.




Figure 2: Scattering diagrams showing the interaction of light waves with a diffraction grating.

The diagrams in Figure 2 illustrate the scattering of light waves from a diffraction grating. Subfigure 2a shows the scattering pattern generated by many grating lines, while Subfigure 2b focuses on the scattering from two grating lines, highlighting the path difference responsible for constructive or destructive interference. These figures provide a foundational understanding of how light interacts with a grating to produce its characteristic diffraction pattern [2].

METHODOLOGY



Figure 3: Experimental setup for the diffraction grating spectrometer. The setup includes the spectrometer (center), a light source (right) emitting the spectrum, and a telescope for observing and measuring diffraction angles. The diffraction grating is positioned between the collimator and the telescope to disperse the light into spectral lines.

The experiment involves using a diffraction grating spectrometer to measure the wavelengths of sodium, mercury, helium, cadmium, and hydrogen. The spectrometer is calibrated by first adjusting the telescope to focus on parallel rays and aligning the collimator for parallel light emission. A diffraction grating is then placed perpendicular to the axis of the collimator on the spectrometer table, ensuring proper alignment. The sodium lamp is switched on, and the telescope crosshairs are set on the central image to establish a zero-angle reference. The grating table is leveled, and the angles of diffraction maxima (θ) on both sides of the light spectrum are measured using the Vernier scale. The measured angles are used in the grating equation:

$$m\lambda = d\sin\theta$$
.

where d is the grating spacing, calculated as

$$d=\frac{1}{N},$$

with N being the number of lines per millimeter on the grating. m is the order of the diffrac-

tion maxima, and λ is the wavelength. The procedure is repeated for mercury, helium, cadmium, and hydrogen discharge lamps. Data analysis includes computing the standard deviation of the wavelengths, comparing experimental wavelengths to standard values, and performing error analysis [3].

DATA ANALYSIS

All data, calculation and programming python code of the experments are attached in the Appendices.

Table 1 presents the spectral line data for various elements, including mercury (Hg), helium (He), cadmium (Cd), hydrogen (H_2) , and sodium (Na).

Element	Color	m	$\theta_L (\deg)$	θ_R (deg)	θ_L (rad)	θ_R (rad)	$\bar{\theta}$ (rad)
Hg	White	0	11°11'30"	11°24'30"	0.195331	0.199113	0.197222
Hg	Purple	1	25°54'30"	26°13'0"	0.452186	0.457567	0.454876
Hg	Green	1	29°47'30"	30°8'30"	0.519963	0.526071	0.523017
Hg	Yellow	1	30°52'50"	31°14'0"	0.538967	0.545125	0.542046
Hg	Purple	2	41°40'50"	42°6'0"	0.727463	0.734784	0.731123
Hg	Green	2	50°39'0"	50°56'30"	0.884009	0.889100	0.886555
Hg	Yellow	2	53°19'0"	53°44'30"	0.930551	0.937969	0.934260
Не	Orange	0	11°9'10"	11°22'40"	0.194653	0.198580	0.196616
Не	Blue	1	26°36'0"	26°54'0"	0.464258	0.469494	0.466876
Не	Cyan	1	28°20'0"	28°38'20"	0.494510	0.499843	0.497176
Не	Yellow	1	31°10'30"	31°27'30"	0.544106	0.549051	0.546579
Cd	Cyan	0	11°17'30"	11°32'0"	0.197077	0.201295	0.199186
Cd	Blue	1	27°17'0"	27°38'10"	0.476184	0.482341	0.479263
Cd	Light Blue	1	27°45'0"	28°4'10"	0.484329	0.489904	0.487117
Cd	Green	1	28°47'10"	29°8'10"	0.502412	0.508521	0.505467
H2	Reddish Pink	0	11°7'20"	11°22'10"	0.194119	0.198434	0.196277
H2	Blue	1	25°52'0"	26°10'40"	0.451458	0.456888	0.454173
H2	Cyan	1	26°40'50"	27°29'40"	0.465664	0.479869	0.472766
H2	Red	1	33°39'20"	34°49'40''	0.587400	0.607859	0.597630
Na	Orange	0	11°8'40"	11°23'10"	0.194507	0.198725	0.196616
Na	Orange	1	31°11'0"	31°25'20"	0.544252	0.548421	0.546337

Table 1: Spectral Line Data for Various Elements with $\bar{\theta}$ in Radians.

The data in **Table** 1 highlight clear trends in the diffraction angles $\bar{\theta}$ across multiple spectral lines for various elements. These angles increase with both the diffraction order (m) and the wavelength of the spectral line, which is consistent with the diffraction grating equation:

$$d\sin\theta = m\lambda$$
.

For mercury (Hg), the yellow spectral line at first-order diffraction (m=1) has a $\bar{\theta}$ value of 0.542046 rad, while at second-order diffraction (m=2), the angle increases to 0.934260 rad. This demonstrates that higher diffraction orders result in larger angular deviations. For helium (He), the angular measurements of the blue, cyan, and yellow lines show a systematic increase in $\bar{\theta}$ with increasing wavelength.

Similar trends are observed for cadmium (Cd), where the first-order blue and green lines exhibit increasing values of $\bar{\theta}$ at 0.479263 rad and 0.505467 rad, respectively. For hydrogen (H_2), the red line has the largest angle among the observed spectral lines, with $\bar{\theta}=0.597630$ rad, reflecting its longer wavelength compared to the blue and cyan lines. The sodium (Na) spectral data further support these observations. The first-order orange line exhibits a $\bar{\theta}$ value of 0.546337 rad, which is consistent with its longer wavelength.

Overall, the results validate the relationship between the diffraction angle and the wavelength of light. The consistency across elements and spectral lines demonstrates the reliability of the experimental setup and highlights the effectiveness of diffraction gratings for precise wavelength determination.

Mercury

Table 2 presents the angular measurements, calculated diffraction angles and wavelengths for mercury's spectral lines.

Color	m	$\theta_L (\deg)$	θ_R (deg)	θ_L (rad)	θ_R (rad)	$\bar{\theta}$ (rad)	θ (rad)	λ (nm)
White	0	11°11'30"	11°24'30"	0.195331	0.199113	0.197222	0.000000	0
Purple	1	25°54'30"	26°13'0"	0.452186	0.457567	0.454877	0.257655	431.4845
Green	1	29°47'30"	30°8'30"	0.519963	0.526071	0.523017	0.325795	542.4718
Yellow	1	30°52'50"	31°14'0"	0.538967	0.545125	0.542046	0.344824	572.3992
Purple	2	41°40'50"	42°6'0"	0.727463	0.734784	0.731124	0.533902	862.7308
Green	2	50°39'0"	50°56'30"	0.883801	0.888333	0.886067	0.688845	1076.3605
Yellow	2	53°19'0"	53°44'30"	0.930109	0.937297	0.933703	0.736481	1137.3867

Table 2: Spectral Line Data for Mercury (Hg).

Based on **Table 2** the angular measurements and corresponding diffraction angles for various spectral lines observed in mercury light. The angles θ_L and θ_R were recorded symmetrically on both sides of the light spectrum. These angles were converted into radians, and the average diffraction angle $\bar{\theta}$ was calculated to determine the diffraction angle θ for each spectral line.

The angles θ provided in **Table 2** are calculated relative to the zeroth-order diffraction peak (m=0), which serves as the reference point for all subsequent diffraction angles. The zeroth-order peak represents the central position where no diffraction occurs, corresponding to $\theta=0$. The angles θ_L and θ_R measured for higher orders (m=1 and m=2) are absolute deviations from this reference point. By averaging these values, the mean diffraction angle $\bar{\theta}$ accurately represents the position of the spectral lines relative to the zeroth order.

The results show that the diffraction angle θ increases with increasing wavelength λ , which is consistent with the diffraction grating equation $d \sin \theta = m\lambda$. For instance, at m = 1, the diffraction angles for the purple, green and yellow lines are 0.257655 rad, 0.325795 rad and 0.344824 rad, corresponding to wavelengths of 431.4845 nm, 542.4718 nm and 572.3992 nm, respectively. For higher orders, such as m = 2, the angles increase further, as seen for the purple, green and yellow lines. This trend validates the theoretical expectation that longer

wavelengths result in greater diffraction angles.

Table 3 compares the wavelengths calculated from the diffraction measurements with the theoretical (real) values of mercury's spectral lines, along with their percentage difference.

Color	m	Calculated Wavelength (nm)	Real Wavelength (nm)	Percentage Difference (%)
Purple	1	431.4845	435.8328	1.00
Green	1	541.9718	546.0735	0.75
Yellow	1	572.3992	567.7105	0.83

Table 3: Mercury spectrum data with calculated and real wavelengths, and their percentage differences.[4][5]

Based on **Table 3**, a comparison between the calculated wavelengths obtained from the diffraction measurements and the theoretical (real) values of mercury's spectral lines, along with their percentage differences. For purple light, the calculated wavelength is **431.4845 nm**, which is within 1.00% of the real value 435.8328 nm. For green light, the calculated wavelength is **541.9718 nm**, with a minor deviation of 0.75% from the theoretical value 546.0735 nm. Similarly, for yellow light, the calculated wavelength of **572.3992 nm** shows a small difference of 0.83% compared to the real value 567.7105 nm.

Helium

Table 4 presents the angular measurements, calculated diffraction angles, and wavelengths for helium's spectral lines.

Color	m	$\theta_L (\deg)$	θ_R (deg)	θ_L (rad)	θ_R (rad)	$\bar{\theta}$ (rad)	θ (rad)	λ (nm)
Orange	0	11°9'10"	11°22'40"	0.194653	0.198580	0.196616	0.0000	0.0000
Blue	1	26°36'0"	26°54'0"	0.464258	0.469494	0.466876	0.2703	452.0885
Cyan	1	28°20'0"	28°38'20"	0.494510	0.499843	0.497176	0.3006	501.3204
Yellow	1	31°10'30"	31°27'30"	0.544106	0.549051	0.546579	0.3500	580.5810

Table 4: Spectral Line Data for Helium (He).

Based on **Table 4**, at m = 1, the diffraction angles for blue, cyan, and yellow light are **0.2703 rad**, **0.3006 rad**, and **0.3500 rad**, corresponding to wavelengths of **452.0885 nm**, **501.3204 nm**, and **580.5810 nm**, respectively. These values confirm that longer wavelengths result in larger diffraction angles, as expected.

Table 5 compares the wavelengths calculated from the diffraction measurements with the theoretical (real) values of helium's spectral lines, along with their percentage differences.

Color	m	Calculated Wavelength (nm)	Real Wavelength (nm)	Percentage Difference (%)
Blue	1	452.0885	447.1479	1.10
Cyan	1	501.3204	501.5678	0.05
Yellow	1	580.5810	587.5615	1.19

Table 5: Helium spectrum data with calculated and real wavelengths, and their percentage differences.[6][7]

Based on **Table 5**, the calculated wavelengths closely match the real (theoretical) values, with small percentage differences observed. For blue light, the calculated wavelength of **452.0885** nm deviates by 1.10% from the real value 447.1479 nm. For cyan light, the calculated wavelength of **501.3204** nm shows excellent agreement with the theoretical value 501.5678 nm, with a negligible percentage difference of 0.05%. Similarly, the yellow light has a calculated wavelength of **580.5810** nm, which deviates by 1.19% from the real value 587.5615 nm.

Cadmium

Table 6 presents the angular measurements, calculated diffraction angles, and wavelengths for cadmium's spectral lines.

Color	m	$\theta_L (\deg)$	θ_R (deg)	θ_L (rad)	θ_R (rad)	$\bar{\theta}$ (rad)	θ (rad)	λ (nm)
Cyan	0	11°17'30"	11°32'0"	0.197077	0.201295	0.199186	0.0000	0.0000
Blue	1	27°17'0"	27°38'10"	0.476184	0.482341	0.479263	0.2801	468.0873
Light Blue	1	27°45'0"	28°4'10"	0.484329	0.489904	0.487117	0.2879	480.8539
Green	1	28°47'10"	29°8'10"	0.502412	0.508521	0.505467	0.3063	510.5651

Table 6: Spectral Line Data for Cadmium (Cd).

Based on **Table** 6, at m = 1, the diffraction angles for blue, light blue, and green spectral lines are **0.2801** rad, **0.2879** rad, and **0.3063** rad, corresponding to wavelengths of **468.0873** nm, **480.8539** nm, and **510.5651** nm, respectively. This observation aligns with the theoretical relationship where the diffraction angle increases with increasing wavelength.

Table 7 compares the wavelengths calculated from the diffraction measurements with the theoretical (real) values of cadmium's spectral lines, along with their percentage differences.

Color	m	Calculated Wavelength (nm)	Real Wavelength (nm)	Percentage Difference (%)
Blue	1	468.0873	467.8149	0.06
Light Blue	1	480.8539	479.9912	0.18
Green	1	510.5651	508.5822	0.39

Table 7: Cadmium spectrum data with calculated and real wavelengths, and their percentage differences.[8][9]

Based on **Table 7**, the calculated wavelengths closely match the real (theoretical) values of cadmium's spectral lines, with small percentage differences observed. For blue light, the calculated wavelength of **468.0873 nm** deviates by only 0.06% from the real value 467.8149 nm. Similarly, for light blue light, the calculated wavelength of **480.8539 nm** has a small difference of 0.18% compared to the real value 479.9912 nm. The green light shows a calculated wavelength of **510.5651 nm**, with a deviation of 0.39% from the real value 508.5822 nm.

Hydrogen

Table 8 presents the angular measurements, calculated diffraction angles, and wavelengths for hydrogen's spectral lines.

Color	m	$\theta_L (\deg)$	θ_R (deg)	θ_L (rad)	θ_R (rad)	$\bar{\theta}$ (rad)	θ (rad)	λ (nm)
Reddish Pink	0	11°7'20"	11°22'10"	0.194119	0.198434	0.196277	0.0000	0.0000
Blue	1	25°52'0"	26°10'40"	0.451458	0.456888	0.454173	0.2579	431.8801
Cyan	1	26°40'50"	27°29'40"	0.465664	0.479869	0.472766	0.2765	462.2460
Red	1	33°39'20"	34°49'40''	0.587400	0.607859	0.597630	0.4014	661.5247

Table 8: Spectral Line Data for Hydrogen (H_2) .

Based on **Table 8**, at m = 1, the diffraction angles for blue, cyan, and red spectral lines are **0.2579** rad, **0.2765** rad, and **0.4014** rad, corresponding to wavelengths of **431.8801** nm, **462.2460** nm, and **661.5247** nm, respectively. These results clearly show that longer wavelengths correspond to larger diffraction angles, which is consistent with theoretical predictions.

Table 9 compares the calculated wavelengths with the theoretical (real) values of hydrogen's spectral lines and their percentage differences.

Color	m	Calculated Wavelength (nm)	Real Wavelength (nm)	Percentage Difference (%)
Blue	1	431.8801	434.0462	0.50
Cyan	1	462.2460	486.1362	4.91
Red	1	661.5247	656.2852	0.80

Table 9: Hydrogen spectrum data with calculated and real wavelengths, and their percentage differences.[10][11]

Based on **Table** 9, the calculated wavelengths are compared with the theoretical values for hydrogen's spectral lines. The results show that the blue light has a calculated wavelength of **431.8801** nm, which deviates by 0.50% from the real value of 434.0462 nm. For cyan light, the calculated wavelength is **462.2460** nm, with a larger percentage difference of 4.91% compared to the theoretical value 486.1362 nm. The red light has a calculated wavelength of **661.5247** nm, showing a small deviation of 0.80% from the real value 656.2852 nm.

Sodium

Table 10 presents the angular measurements, calculated diffraction angles, and wavelengths for sodium's spectral lines.

Color	m	θ_L (deg)	θ_R (deg)	θ_L (rad)	θ_R (rad)	$\bar{\theta}$ (rad)	θ (rad)	λ (nm)
Orange	0	11°8'40"	11°23'10"	0.194507	0.198725	0.196616	0.00000	0.0000
Orange	1	31°11'0"	31°25'20"	0.544252	0.548421	0.546337	0.34972	580.1954

Table 10: Spectral Line Data for Sodium (Na).

Based on **Table 10**, at m = 1, the diffraction angle for the orange line is **0.34972 rad**, corresponding to a calculated wavelength of **580.1954 nm**. This observation confirms the relationship where longer wavelengths are associated with larger diffraction angles.

Table 11 compares the calculated wavelength with the theoretical (real) value of sodium's spectral line and the percentage difference.

Color	m	Calculated Wavelength (nm)	Real Wavelength (nm)	Percentage Difference (%)
Orange	1	580.1954	588.995	1.49

Table 11: Sodium spectrum data with calculated and real wavelengths, and their percentage differences.[12][13]

Based on **Table 11**, the calculated wavelength for the orange line is **580.1954 nm**, which deviates by 1.49% from the real wavelength of 588.995 nm. This small discrepancy can be attributed to minor uncertainties in the angular measurements, instrumental resolution, or imperfections in the diffraction grating. Despite the minor difference, the experimental results closely align with the theoretical value, validating the diffraction grating method as an accurate tool for measuring spectral line wavelengths.

DISCUSSION

In this experiment, a diffraction grating spectrometer was employed to measure the wavelengths of light emitted by sodium, mercury, helium, cadmium, and hydrogen. The grating equation, $\lambda = m\lambda$, was utilized to calculate the wavelengths, where θ is the diffraction angle, d is the grating spacing, and m represents the diffraction order. This equation provided the foundation for analyzing the diffraction patterns observed.

In the case of **mercury** (Hg), the first-order spectral lines for purple, green, and yellow light were measured. The calculated wavelengths were 431.4845 nm, 541.9718 nm, and 572.3992 nm, compared to their respective theoretical values of 435.8328 nm, 546.0735 nm, and 567.7105 nm. The percentage differences were 1.00%, 0.75%, and 0.83%, respectively. These results confirm the experiment's precision, with minor discrepancies arising from small angular reading errors.

For **helium** (*He*), the blue, cyan, and yellow spectral lines were analyzed. The calculated wavelengths were 452.0885 nm, 501.3204 nm, and 580.5810 nm, compared to their theoretical values of 447.1479 nm, 501.5678 nm, and 587.5615 nm. The percentage differences were 1.10%, 0.05%, and 1.19%, respectively. The cyan line showed an exceptional agreement with the theoretical value, confirming the diffraction grating's reliability, while minor deviations for the other lines reflected systematic measurement uncertainties.

For **cadmium** (Cd), the blue, light blue, and green lines were observed at m=1. The calculated wavelengths were $468.0873 \,\mathrm{nm}$, $480.8539 \,\mathrm{nm}$, and $510.5651 \,\mathrm{nm}$, compared to their theoretical values of $467.8149 \,\mathrm{nm}$, $479.9912 \,\mathrm{nm}$, and $508.5822 \,\mathrm{nm}$. The corresponding percentage differences were 0.06%, 0.18%, and 0.39%, indicating high precision in the measurements.

For **hydrogen** (H_2), the blue, cyan, and red spectral lines were analyzed. The calculated wavelengths were 431.8801 nm, 462.2460 nm, and 661.5247 nm, with theoretical values of 434.0462 nm, 486.1362 nm, and 656.2852 nm, respectively. The percentage differences were 0.50% for the blue line, 4.91% for the cyan line, and 0.80% for the red line. The larger deviation for the cyan line likely resulted from difficulties in observing faint spectral lines and alignment errors.

For **sodium** (Na), the measured orange spectral line at m=1 produced a calculated wavelength of 580.1954 nm, compared to the theoretical value of 588.995 nm, yielding a percentage difference of 1.49%. This small deviation reflects the method's accuracy, with slight errors attributed to alignment issues and instrumental resolution.

Various sources of error were identified during the experiment. Instrumental factors included potential imperfections or damage to the diffraction grating, which could affect the uniformity of line spacing and introduce inaccuracies in measured angles. The presence of dust on the diffraction grating's lens may have caused unwanted reflections of light, interfering with the clarity and precision of the diffraction pattern. Misalignment of the collimator and telescope also posed a risk to angle measurements, impacting the accuracy of the calculated wavelengths. The inability to adjust the focus of the spectrometer to the actual position of the spectral lines could have further affected the accuracy of the angular readings, leading to minor discrepancies between the calculated and theoretical wavelengths. The low intensity of the light sources contributed to the dimness of the observed spectral lines, making it difficult to clearly identify the boundaries of the light spectrum. This limitation could have introduced errors in measuring the diffraction angles. Environmental influences, such as temperature fluctuations and air currents, altered the refractive index, potentially affecting the diffraction patterns. Human limitations, particularly the eye's reduced sensitivity to dim spectral lines in the red and violet regions, further complicated precise observations.

To address these challenges, several improvements were adopted. Cleaning the dust on diffraction grating glass by using alcohol and cotton swab. Regular calibration of the spectrometer using light sources with known wavelengths ensured systematic errors were identified and corrected. Experiments were conducted in controlled environments to minimize external disturbances like temperature variations and air currents. Sensitive detection equipment, such as photomultiplier tubes or CCDs (Charge-coupled Device), enhanced the visibility of dim spectral lines, compensating for human limitations. Increasing the number of slits in the diffraction grating was another effective strategy, as it sharpened the diffraction maxima and improved the resolution and precision of wavelength measurements.

CONCLUSION

For the Mercury spectrum, the wavelengths are 431.48 nm for purple, 541.97 nm for green, and 572.40 nm for yellow.

For the Helium spectrum, the wavelengths are $452.09\,\mathrm{nm}$ for blue, $501.32\,\mathrm{nm}$ for cyan, and $580.58\,\mathrm{nm}$ for yellow.

For the Cadmium spectrum, the wavelengths are 468.09 nm for blue, 480.85 nm for light blue, and 510.57 nm for green.

For the Hydrogen spectrum, the wavelengths are 431.88 nm for blue, 462.25 nm for cyan, and 661.52 nm for red.

For the Sodium spectrum, the wavelengths are 580.20 nm for orange.

REFERENCES

- [1] P. Scientific. Educational Spectrophotometer System Manual. https://cdn.pasco.com/product_document/Educational-Spectrophotometer-System-Manual-OS-8539.pdf. Accessed: November 3, 2024. 2017.
- [2] S. Cortland. *Physics Laboratory Manual: Experiment 1*. https://web.cortland.edu/wheeler/phy150/Laboratory/Laboratory1.pdf. Accessed: November 3, 2024. Nov. 2024.
- [3] Universiti Sains Malaysia (USM). 20S2 Diffraction Grating Spectrometer Manual. [Accessed: Dec. 8, 2024]. Dec. 2024.
- [4] U. contributions. Second order spectra vs. first order. https://physics.stackexchange.com/questions/284394/second-order-spectra-vs-first-order. Accessed: 2024-12-15.
- [5] N. I. of Standards and T. (NIST). *Atomic Spectra Database Mercury Wavelengths*. https://physics.nist.gov/PhysRefData/Handbook/Tables/mercurytable2.htm. Accessed: 2024-12-15.
- [6] S. P. Library. Helium emission and absorption spectra. https://www.sciencephoto.com/media/673903/view/helium-emission-and-absorption-spectra. Accessed: 2024-12-15.
- [7] N. I. of Standards and T. (NIST). *Atomic Spectra Database Helium Wavelengths*. https://physics.nist.gov/PhysRefData/Handbook/Tables/heliumtable2.htm. Accessed: 2024-12-15.
- [8] S. P. Library. *Emission Spectrum of Cadmium*. https://www.sciencephoto.com/media/1674/view/emission-spectrum-of-cadmium. Accessed: 2024-12-15.
- [9] N. I. of Standards and T. (NIST). Atomic Spectra Database Cadmium Wavelengths. https://physics.nist.gov/PhysRefData/Handbook/Tables/cadmiumtable2_a.htm. Accessed: 2024-12-15.
- [10] U. of Texas at Austin. *Hydrogen Atom Line Spectra*. https://ch301.cm.utexas.edu/section2.php?target=atomic/H-atom/line-spectra.html. Accessed: 2024-12-15.
- [11] N. I. of Standards and T. (NIST). *Atomic Spectra Database Hydrogen Wavelengths*. https://physics.nist.gov/PhysRefData/Handbook/Tables/hydrogentable2.htm. Accessed: 2024-12-15.

- [12] L. Gallery. Sodium Vapor Emission Spectrum. https://www.lighting-gallery.net/gallery/displayimage.php?album=3248&pos=46&pid=110444. Accessed: 2024-12-15.
- [13] N. I. of Standards and T. (NIST). *Atomic Spectra Database Sodium Wavelengths*. https://physics.nist.gov/PhysRefData/Handbook/Tables/sodiumtable2.htm. Accessed: 2024-12-15.

APPENDICES

Data

```
import math
    import pandas as pd
    # Function to convert degrees, minutes, and seconds to radians
    def dms_to_radians(degrees, minutes, seconds):
        total_degrees = degrees + minutes / 60 + seconds / 3600
        return math.radians(total_degrees)
    # Data for the table
    data = [
        ["Hg", "White", 0, (11, 11, 30), (11, 24, 30)],
        ["Hg", "Purple", 1, (25, 54, 30), (26, 13, 0)],
        ["Hg", "Green", 1, (29, 47, 30), (30, 8, 30)],
        ["Hg", "Yellow", 1, (30, 52, 50), (31, 14, 0)],
        ["Hg", "Purple", 2, (41, 40, 50), (42, 6, 0)],
        ["Hg", "Green", 2, (50, 39, 0), (50, 56, 30)],
        ["Hg", "Yellow", 2, (53, 19, 0), (53, 44, 30)],
        ["He", "Orange", 0, (11, 9, 10), (11, 22, 40)],
        ["He", "Blue", 1, (26, 36, 0), (26, 54, 0)],
        ["He", "Cyan", 1, (28, 20, 0), (28, 38, 20)],
        ["He", "Yellow", 1, (31, 10, 30), (31, 27, 30)],
21
        ["Cd", "Cyan", 0, (11, 17, 30), (11, 32, 0)],
        ["Cd", "Blue", 1, (27, 17, 0), (27, 38, 10)],
        ["Cd", "Light Blue", 1, (27, 45, 0), (28, 4, 10)],
        ["Cd", "Green", 1, (28, 47, 10), (29, 8, 10)],
        ["H2", "Reddish Pink", 0, (11, 7, 20), (11, 22, 10)],
        ["H2", "Blue", 1, (25, 52, 0), (26, 10, 40)],
        ["H2", "Cyan", 1, (26, 40, 50), (27, 29, 40)],
        ["H2", "Red", 1, (33, 39, 20), (34, 49, 40)],
        ["Na", "Orange", 0, (11, 8, 40), (11, 23, 10)],
        ["Na", "Orange", 1, (31, 11, 0), (31, 25, 20)],
32
33
    # Process the data
34
    processed_data = []
35
36
    for element, color, m, theta_L, theta_R in data:
37
        theta_L_rad = dms_to_radians(*theta_L)
        theta_R_rad = dms_to_radians(*theta_R)
```

Output

	P3		0.3	TT 1.		+- D (d) This T	(4)
1		ment		The	_	ta_R (deg) Theta_L	
2	0	Hg	White			(11, 24, 30)	0.195331
3	1	Hg	Purple	1		(26, 13, 0)	0.452186
4	2	Hg	Green	1	(29, 47, 30)	(30, 8, 30)	0.519963
5	3	Hg 	Yellow	1	(30, 52, 50)	(31, 14, 0)	0.538967
6	4	Hg	Purple	2	(41, 40, 50)	(42, 6, 0)	0.727463
7	5	Hg 	Green	2	(50, 39, 0)	(50, 56, 30)	0.884009
8	6	Hg	Yellow	2	(53, 19, 0)	(53, 44, 30)	0.930551
9	7	He	Orange	0	(11, 9, 10)	(11, 22, 40)	0.194653
10	8	He	Blue	1	(26, 36, 0)	(26, 54, 0)	0.464258
11	9	He	Cyan	1	(28, 20, 0)	(28, 38, 20)	0.494510
12	10	He	Yellow	1	(31, 10, 30)	(31, 27, 30)	0.544106
13	11	Cd	Cyan	0	(11, 17, 30)	(11, 32, 0)	0.197077
14	12	Cd	Blue	1	(27, 17, 0)	(27, 38, 10)	0.476184
15	13	Cd	Light Blue	1		(28, 4, 10)	0.484329
16	14	Cd	Green	1		(29, 8, 10)	0.502412
17	15	Н2	Reddish Pink	0	(11, 7, 20)	(11, 22, 10)	0.194119
18	16	Н2	Blue	1	(25, 52, 0)	(26, 10, 40)	0.451458
19	17	Н2	Cyan	1	(26, 40, 50)	(27, 29, 40)	0.465664
20	18	Н2	Red	1	(33, 39, 20)	(34, 49, 40)	0.587400
21	19	Na	Orange	0	(11, 8, 40)	(11, 23, 10)	0.194507
22	20	Na	Orange	1	(31, 11, 0)	(31, 25, 20)	0.544252
23							
24		Theta_R	(rad) Delta	The	eta (rad)		
25	0		199113		0.197222		
26	1	0.	457567		0.454876		
27	2		526071		0.523017		
28	3	0.	545125		0.542046		
29	4	0.	734784		0.731123		
30	5		889100		0.886555		
31	6	0.	937969		0.934260		
32	7	0.	198580		0.196616		
33	8	0.	469494		0.466876		
34	9	0.	499843		0.497176		
35	10	0.	549051		0.546579		
36	11	0.	201295		0.199186		
37	12	0.	482341		0.479263		
38	13	0.	489904		0.487117		
39	14	0.	508521		0.505467		
40	15	0.	198434		0.196277		

41	16	0.456888	0.454173	
42	17	0.479869	0.472766	
43	18	0.607859	0.597630	
44	19	0.198725	0.196616	
45	20	0.548421	0.546337	

Mercury

```
import math
    import pandas as pd
    # Data for the table
    data = [
        ["White", 0, (11, 11, 30), (11, 24, 30), 0.195331, 0.199113, 0.197222],
        ["Purple", 1, (25, 54, 30), (26, 13, 0), 0.452186, 0.457567, 0.454877],
        ["Green", 1, (29, 47, 30), (30, 8, 30), 0.519963, 0.526071, 0.523017],
        ["Yellow", 1, (30, 52, 50), (31, 14, 0), 0.538967, 0.545125, 0.542046],
        ["Purple", 2, (41, 40, 50), (42, 6, 0), 0.727463, 0.734784, 0.731124],
        ["Green", 2, (50, 39, 0), (50, 56, 30), 0.883801, 0.888333, 0.886067],
        ["Yellow", 2, (53, 19, 0), (53, 44, 30), 0.930109, 0.937297, 0.933703],
13
   ]
14
    # Create a DataFrame
    columns = [
16
        "Color", "m", "Theta_L (deg)", "Theta_R (deg)", "Theta_L (rad)", "Theta_R
      (rad)", "Delta Theta (rad)"
    df = pd.DataFrame(data, columns=columns)
19
20
    # Extract Delta Theta for m = 0
21
   delta_theta_m0 = df[df["m"] == 0]["Delta Theta (rad)"].iloc[0]
23
    # Calculate theta for each color
    df["Theta (rad)"] = df["Delta Theta (rad)"] - delta_theta_m0
26
    # Grating spacing in nanometers
27
    d_{inch} = 1 / 15000 # d in inches
28
    d_meter = d_inch * 0.0254 # Convert d to meters
    d_nm = d_meter * 1e9 # Convert d to nanometers
30
    # Calculate d * sin(Theta) in nanometers
    df["d sin(Theta) (nm)"] = (d_nm * df["Theta (rad)"].apply(math.sin)).round(4)
    # Display the updated DataFrame
    print(df)
```

Output

```
White 0 (11, 11, 30)
                            (11, 24, 30)
                                              0.195331
                                                            0.199113
   1 Purple 1 (25, 54, 30)
                            (26, 13, 0)
                                              0.452186
                                                            0.457567
     Green 1 (29, 47, 30)
                            (30, 8, 30)
                                                            0.526071
                                              0.519963
   3 Yellow 1 (30, 52, 50)
                            (31, 14, 0)
                                              0.538967
                                                            0.545125
                             (42, 6, 0)
   4 Purple 2 (41, 40, 50)
                                              0.727463
                                                            0.734784
                (50, 39, 0)
                            (50, 56, 30)
                                              0.883801
                                                            0.888333
   5
      Green 2
     Yellow 2
                 (53, 19, 0)
                             (53, 44, 30)
                                              0.930109
                                                            0.937297
      Delta Theta (rad) Theta (rad) d sin(Theta) (nm)
10
11
   0
              0.197222
                          0.00000
                                             0.0000
              0.454877
                          0.257655
                                           431.4845
12
   1
              0.523017
                          0.325795
                                           541.9718
13
   2
   3
              0.542046
                          0.344824
                                           572.3992
14
   4
              0.731124
                          0.533902
                                           861.7308
15
              0.886067
                          0.688845
                                          1076.3605
16
   5
                                          1137.3867
              0.933703
                          0.736481
```

Mercury percentage difference

```
import pandas as pd
    # Provided data for m = 1
    data = {
       "Color": ["Purple", "Green", "Yellow"],
       "m": [1, 1, 1],
       "Calculated Wavelength (nm)": [431.4845, 541.9718, 572.3992],
   }
    \# Add the real wavelength values provided by the user for m=1
   real_values = [435.8328, 546.0735, 567.7105] # Replace with your actual
     values
   data["Real Wavelength (nm)"] = real_values
    # Create a DataFrame
15
   df = pd.DataFrame(data)
    # Calculate the percentage difference
    df["Percentage Difference (%)"] = (
        abs(df["Calculated Wavelength (nm)"] - df["Real Wavelength (nm)"]) / df["
    Real Wavelength (nm)"] * 100
   ).round(2)
21
    # Show only the required columns
   result_df = df[["Color", "m", "Calculated Wavelength (nm)", "Real Wavelength
     (nm)", "Percentage Difference (%)"]]
    # Display the result DataFrame
   print(result_df)
```

Output

```
Color m Calculated Wavelength (nm) Real Wavelength (nm) \
0 Purple 1
                             431.4845
                                                435.8328
                            541.9718
                                               546.0735
1 Green 1
                            572.3992
                                               567.7105
2 Yellow 1
 Percentage Difference (%)
                      1.00
0
                      0.75
1
                      0.83
```

Helium

```
import math
    import pandas as pd
    # Function to convert degrees, minutes, and seconds to radians
    def dms_to_radians(degrees, minutes, seconds):
        total_degrees = degrees + minutes / 60 + seconds / 3600
        return math.radians(total_degrees)
    # Data for the table
    data = [
        ["Orange", 0, (11, 9, 10), (11, 22, 40)],
        ["Blue", 1, (26, 36, 0), (26, 54, 0)],
        ["Cyan", 1, (28, 20, 0), (28, 38, 20)],
        ["Yellow", 1, (31, 10, 30), (31, 27, 30)],
15
16
    # Process the data
    processed_data = []
    for color, m, theta_L_dms, theta_R_dms in data:
        theta_L_rad = dms_to_radians(*theta_L_dms)
        theta_R_rad = dms_to_radians(*theta_R_dms)
        delta_theta_rad = (theta_L_rad + theta_R_rad) / 2
        {\tt processed\_data.append([color, m, theta\_L\_dms, theta\_R\_dms, theta\_L\_rad,}
     theta_R_rad, delta_theta_rad])
    # Create a DataFrame
    columns = [
        "Color", "m", "Theta_L (deg)", "Theta_R (deg)", "Theta_L (rad)", "Theta_R
      (rad)", "Delta Theta (rad)"
    df = pd.DataFrame(processed_data, columns=columns)
29
    # Extract Delta Theta for m = 0
    delta_theta_m0 = df[df["m"] == 0]["Delta Theta (rad)"].iloc[0]
    # Calculate theta for each color
   df["Theta (rad)"] = df["Delta Theta (rad)"] - delta_theta_m0
35
36
37
    # Grating spacing in nanometers
   d_{inch} = 1 / 15000 # d in inches
```

```
d_meter = d_inch * 0.0254 # Convert d to meters

d_nm = d_meter * 1e9 # Convert d to nanometers

# Calculate d * sin(Theta) in nanometers

df["d sin(Theta) (nm)"] = (d_nm * df["Theta (rad)"].apply(math.sin)).round(4)

# Display the updated DataFrame

print(df)
```

```
0 Orange 0 (11, 9, 10) (11, 22, 40)
                                         0.194653
                                                     0.198580
      Blue 1 (26, 36, 0) (26, 54, 0)
                                         0.464258
                                                     0.469494
  1
       Cyan 1 (28, 20, 0) (28, 38, 20)
                                                     0.499843
                                         0.494510
   3 Yellow 1 (31, 10, 30) (31, 27, 30)
                                         0.544106
                                                     0.549051
     Delta Theta (rad) Theta (rad) d sin(Theta) (nm)
            0.196616
                       0.000000
                                        0.0000
   1
            0.466876
                       0.270259
                                      452.0885
10
  2
            0.497176
                       0.300560
                                      501.3204
            0.546579
                       0.349963
                                      580.5810
   3
```

Helium percentage difference

```
import pandas as pd
    \# Provided data for m = 1
    data = {
        "Color": ["Blue", "Cyan", "Yellow"],
        "m": [1, 1, 1],
        "Calculated Wavelength (nm)": [452.0885,501.3204,580.5810]
   }
   # Real wavelength values for m = 1
   real_values = [447.1479,501.5678,587.5615] # Replace with actual values
    # Create a DataFrame
13
   df = pd.DataFrame(data)
15
    # Add real wavelength values
16
    df["Real Wavelength (nm)"] = real_values
19
    # Calculate the percentage difference
   df["Percentage Difference (%)"] = (
20
        abs(df["Calculated Wavelength (nm)"] - df["Real Wavelength (nm)"]) / df["
     Real Wavelength (nm)"] * 100
   ).round(2)
    # Display the DataFrame
   print(df)
```

```
import pandas as pd
    \# Provided data for m = 1
    data = {
        "Color": ["Blue", "Cyan", "Yellow"],
        "m": [1, 1, 1],
        "Calculated Wavelength (nm)": [452.0885,501.3204,580.5810]
   }
   # Real wavelength values for m = 1
10
   real_values = [447.1479,501.5678,587.5615] # Replace with actual values
    # Create a DataFrame
13
   df = pd.DataFrame(data)
15
    # Add real wavelength values
16
17
    df["Real Wavelength (nm)"] = real_values
19
    # Calculate the percentage difference
   df["Percentage Difference (%)"] = (
20
        abs(df["Calculated Wavelength (nm)"] - df["Real Wavelength (nm)"]) / df["
     Real Wavelength (nm)"] * 100
   ).round(2)
    # Display the DataFrame
   print(df)
```

Cadmium

```
import math
    import pandas as pd
    # Function to convert degrees, minutes, and seconds to radians
    def dms_to_radians(degrees, minutes, seconds):
        total_degrees = degrees + minutes / 60 + seconds / 3600
        return math.radians(total_degrees)
    # Data for the table
    data = [
        ["Cyan", 0, (11, 17, 30), (11, 32, 0)],
        ["Blue", 1, (27, 17, 0), (27, 38, 10)],
        ["Light Blue", 1, (27, 45, 0), (28, 4, 10)],
        ["Green", 1, (28, 47, 10), (29, 8, 10)],
15
16
    # Process the data
    processed_data = []
    for color, m, theta_L_dms, theta_R_dms in data:
        theta_L_rad = dms_to_radians(*theta_L_dms)
        theta_R_rad = dms_to_radians(*theta_R_dms)
        delta_theta_rad = (theta_L_rad + theta_R_rad) / 2
        {\tt processed\_data.append([color, m, theta\_L\_dms, theta\_R\_dms, theta\_L\_rad,}
     theta_R_rad, delta_theta_rad])
    # Create a DataFrame
    columns = [
        "Color", "m", "Theta_L (deg)", "Theta_R (deg)", "Theta_L (rad)", "Theta_R
      (rad)", "Delta Theta (rad)"
28
    df = pd.DataFrame(processed_data, columns=columns)
29
    # Extract Delta Theta for m = 0
    delta_theta_m0 = df[df["m"] == 0]["Delta Theta (rad)"].iloc[0]
    # Calculate theta for each color
   df["Theta (rad)"] = df["Delta Theta (rad)"] - delta_theta_m0
35
36
    # Grating spacing in nanometers
   d_{inch} = 1 / 15000 # d in inches
```

```
d_meter = d_inch * 0.0254 # Convert d to meters
d_nm = d_meter * 1e9 # Convert d to nanometers

# Calculate d * sin(Theta) in nanometers

df["d sin(Theta) (nm)"] = (d_nm * df["Theta (rad)"].apply(math.sin)).round(4)

# Display the updated DataFrame
print(df)
```

```
Color m Theta_L (deg) Theta_R (deg) Theta_L (rad) Theta_R (rad) \
            Cyan 0 (11, 17, 30) (11, 32, 0)
                                                    0.197077
                                                                  0.201295
            Blue 1 (27, 17, 0) (27, 38, 10)
                                                    0.476184
                                                                  0.482341
   1
   2 Light Blue 1 (27, 45, 0) (28, 4, 10)
                                                    0.484329
                                                                  0.489904
           Green 1 (28, 47, 10) (29, 8, 10)
                                                    0.502412
                                                                  0.508521
      Delta Theta (rad) Theta (rad) d sin(Theta) (nm)
              0.199186
                           0.000000
                                               0.0000
                           0.280077
   1
              0.479263
                                             468.0873
10
   2
              0.487117
                           0.287931
                                             480.8539
              0.505467
                           0.306281
                                             510.5651
   3
```

Cadmium percentage difference

```
import pandas as pd
    # Provided data for m = 1
    data = {
        "Color": ["Blue", "Light Blue", "Green"],
        "m": [1, 1, 1],
        "Calculated Wavelength (nm)": [468.0873, 480.8539, 510.5651]
   }
    # Real wavelength values for m = 1
   real_values = [467.8149, 479.9912, 508.5822] # Replace with actual real
     values
12
    # Create a DataFrame
13
    df = pd.DataFrame(data)
15
    # Add real wavelength values
16
    df["Real Wavelength (nm)"] = real_values
    # Calculate the percentage difference
19
    df["Percentage Difference (%)"] = (
20
        abs(df["Calculated Wavelength (nm)"] - df["Real Wavelength (nm)"]) / df["
21
     Real Wavelength (nm)"] * 100
   ).round(2)
    # Display the result DataFrame
    print(df)
```

```
Color m Calculated Wavelength (nm) Real Wavelength (nm) \setminus
                                     468.0873
                                                            467.8149
         Blue 1
                                     480.8539
1 Light Blue 1
                                                            479.9912
        Green 1
                                     510.5651
                                                            508.5822
   Percentage Difference (%)
                         0.06
0
                         0.18
1
2
                         0.39
```

Hydrogen

```
import math
    import pandas as pd
    # Function to convert degrees, minutes, and seconds to radians
    def dms_to_radians(degrees, minutes, seconds):
        total_degrees = degrees + minutes / 60 + seconds / 3600
        return math.radians(total_degrees)
    # Data for the table
    data = [
        ["Reddish Pink", 0, (11, 7, 20), (11, 22, 10)],
        ["Blue", 1, (25, 52, 0), (26, 10, 40)],
        ["Cyan", 1, (26, 40, 50), (27, 29, 40)],
        ["Red", 1, (33, 39, 20), (34, 49, 40)],
15
16
    # Process the data
    processed_data = []
    for color, m, theta_L_dms, theta_R_dms in data:
        theta_L_rad = dms_to_radians(*theta_L_dms)
        theta_R_rad = dms_to_radians(*theta_R_dms)
        delta_theta_rad = (theta_L_rad + theta_R_rad) / 2
        processed_data.append([color, m, theta_L_dms, theta_R_dms, theta_L_rad,
     theta_R_rad, delta_theta_rad])
    # Create a DataFrame
    columns = [
        "Color", "m", "Theta_L (deg)", "Theta_R (deg)", "Theta_L (rad)", "Theta_R
      (rad)", "Delta Theta (rad)"
    df = pd.DataFrame(processed_data, columns=columns)
29
    # Extract Delta Theta for m = 0
    delta_theta_m0 = df[df["m"] == 0]["Delta Theta (rad)"].iloc[0]
    # Calculate theta for each color
35
    df["Theta (rad)"] = df["Delta Theta (rad)"] - delta_theta_m0
36
    # Grating spacing in nanometers
   d_{inch} = 1 / 15000 # d in inches
```

```
d_meter = d_inch * 0.0254 # Convert d to meters

d_nm = d_meter * 1e9 # Convert d to nanometers

# Calculate d * sin(Theta) in nanometers

df["d sin(Theta) (nm)"] = (d_nm * df["Theta (rad)"].apply(math.sin)).round(4)

# Display the updated DataFrame

print(df)
```

```
Color m Theta_L (deg) Theta_R (deg) Theta_L (rad) Theta_R (rad) \
0 Reddish Pink 0 (11, 7, 20) (11, 22, 10)
                                                   0.194119
                                                                 0.198434
          Blue 1 (25, 52, 0) (26, 10, 40)
                                                   0.451458
                                                                 0.456888
          Cyan 1 (26, 40, 50) (27, 29, 40)
                                                   0.465664
                                                                 0.479869
           Red 1 (33, 39, 20) (34, 49, 40)
                                                   0.587400
                                                                 0.607859
  Delta Theta (rad) Theta (rad) d sin(Theta) (nm)
           0.196277
                        0.000000
                                            0.0000
           0.454173
                        0.257897
                                          431.8801
1
2
           0.472766
                        0.276489
                                          462.2460
           0.597630
                        0.401353
                                          661.5247
3
```

Hydrogen percentage difference

```
import pandas as pd
    \# Provided data for m = 1
    data = {
        "Color": ["Blue", "Cyan", "Red"],
        "m": [1, 1, 1],
        "Calculated Wavelength (nm)": [431.8801, 462.2460, 661.5247]
   }
10
    # Real wavelength values for m = 1
   real_values = [434.0462, 486.1362, 656.2852] # Replace with actual real
     values
    # Create a DataFrame
    df = pd.DataFrame(data)
15
    # Add real wavelength values
16
    df["Real Wavelength (nm)"] = real_values
    # Calculate the percentage difference
19
    df["Percentage Difference (%)"] = (
20
        abs(df["Calculated Wavelength (nm)"] - df["Real Wavelength (nm)"]) / df["
21
     Real Wavelength (nm)"] * 100
   ).round(2)
    # Display the result DataFrame
   print(df)
```

```
Color m Calculated Wavelength (nm) Real Wavelength (nm) \
0 Blue 1
                             431.8801
                                                  434.0462
                             462.2460
1 Cyan 1
                                                  486.1362
  Red 1
                             661.5247
                                                   656.2852
  Percentage Difference (%)
                       0.50
0
                       4.91
1
2
                       0.80
```

Sodium

```
import math
    import pandas as pd
    # Function to convert degrees, minutes, and seconds to radians
    def dms_to_radians(degrees, minutes, seconds):
        total_degrees = degrees + minutes / 60 + seconds / 3600
        return math.radians(total_degrees)
    # Data for the table
    data = [
10
        ["Orange", 0, (11, 8, 40), (11, 23, 10)],
        ["Orange", 1, (31, 11, 0), (31, 25, 20)],
13
14
    # Process the data
    processed_data = []
16
    for color, m, theta_L_dms, theta_R_dms in data:
        theta_L_rad = dms_to_radians(*theta_L_dms)
        theta_R_rad = dms_to_radians(*theta_R_dms)
19
        delta_theta_rad = (theta_L_rad + theta_R_rad) / 2
        {\tt processed\_data.append([color, m, theta\_L\_dms, theta\_R\_dms, theta\_L\_rad,}
     theta_R_rad, delta_theta_rad])
    # Create a DataFrame
    columns = [
        "Color", "m", "Theta_L (deg)", "Theta_R (deg)", "Theta_L (rad)", "Theta_R
      (rad)", "Delta Theta (rad)"
27
    df = pd.DataFrame(processed_data, columns=columns)
28
    # Extract Delta Theta for m = 0
    delta_theta_m0 = df[df["m"] == 0]["Delta Theta (rad)"].iloc[0]
    # Calculate theta for each color
32
    df["Theta (rad)"] = df["Delta Theta (rad)"] - delta_theta_m0
33
35
    # Grating spacing in nanometers
    d_{inch} = 1 / 15000 # d in inches
36
    d_{meter} = d_{inch} * 0.0254 # Convert d to meters
    d_nm = d_meter * 1e9 # Convert d to nanometers
```

```
# Calculate d * sin(Theta) in nanometers

df ["d sin(Theta) (nm)"] = (d_nm * df ["Theta (rad)"].apply(math.sin)).round(4)

# Display the updated DataFrame

print(df)
```

```
Color m Theta_L (deg) Theta_R (deg) Theta_L (rad) Theta_R (rad) \
0 Orange 0 (11, 8, 40) (11, 23, 10) 0.194507 0.198725
1 Orange 1 (31, 11, 0) (31, 25, 20) 0.544252 0.548421

Delta Theta (rad) Theta (rad) d sin(Theta) (nm)
0 0.196616 0.00000 0.0000
7 0.546337 0.34972 580.1954
```

Sodium percentage difference

```
import pandas as pd
   \# Provided data for m = 1
   data = {
       "Color": ["Orange"],
       "m": [1],
       "Calculated Wavelength (nm)": [580.1954]
   }
   \# Real wavelength value for m = 1
   real_values = [588.9950] # Replace with actual real value
   # Create a DataFrame
   df = pd.DataFrame(data)
15
   # Add real wavelength value
16
   df["Real Wavelength (nm)"] = real_values
19
   # Calculate the percentage difference
   df["Percentage Difference (%)"] = (
20
       abs(df["Calculated Wavelength (nm)"] - df["Real Wavelength (nm)"]) / df["
     Real Wavelength (nm)"] * 100
   ).round(2)
   # Display the result DataFrame
   print(df)
```

```
Color m Calculated Wavelength (nm) Real Wavelength (nm) \
0 Orange 1 580.1954 588.995

Percentage Difference (%)
5 O 1.49
```